



Product Datasheet

BEAN1 antibody (orb140147)

Description	Rabbit Polyclonal to BEAN1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	BEAN1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: BEAN. Gene-IDs: 146227. Accessions: NP_001129578.1, Q3B7T3, NP_001171491.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	28626
Uniprot ID	Q3B7T3
NCBI	NP_001184154.1 , NP_001171491.1 , NP_001129578.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...