



Product Datasheet

TGFBRAP1 antibody (orb140153)

Description	Rabbit Polyclonal to TGFBRAP1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	TGFBRAP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: TGFBRAP1. Gene-IDs: 9392. Accessions: Q8WUH2, NP_004248.2, NP_001136093.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	97158
Uniprot ID	Q8WUH2
NCBI	NP_001136093.1 , NP_004248.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...