




Product Datasheet

TRIM38 antibody (orb140156)

Description	Rabbit Polyclonal to TRIM38.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	TRIM38
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: TRIM38. Gene-IDs: 10475. Accessions: NP_006346.1, O00635.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53416
Uniprot ID	O00635
NCBI	NP_006346.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of HeLa cells
(...)