

Product Datasheet

SLC2A9 antibody (orb140163)

Description	Rabbit Polyclonal to SLC2A9.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SLC2A9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: SLC2A9. Gene-IDs: 56606. Accessions: NP_064425.2, NP_001001290.1, Q9NRM0.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	58702
Uniprot ID	Q9NRM0
NCBI	NP_064425.2 , NP_001001290.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

