



Product Datasheet

POU5F2 antibody (orb140165)

Description	Rabbit Polyclonal to POU5F2.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	POU5F2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: POU5F2. Gene-IDs: 134187. Accessions: NP_694948.1, Q8N7G0.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36051
Uniprot ID	Q8N7G0
NCBI	NP_694948.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cell line

...