

Product Datasheet

Doublecortin antibody (orb140166)

Description	Rabbit polyclonal antibody to Doublecortin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DCDC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: DCDC2. Gene-IDs: 51473. Accessions: NP_001182539.1, Q9UHG0, NP_057440.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52834
Uniprot ID	Q9UHG0
NCBI	NP_057440.2 , NP_001182539.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

