

Product Datasheet

MGLL antibody (orb140168)

Description	Rabbit Polyclonal to MGLL.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MGLL
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: MGLL. Gene-IDs: 11343. Accessions: Q99685, NP_001003794.1, NP_009214.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33261
Uniprot ID	Q99685
NCBI	NP_001003794.1 , NP_009214.1 , NP_001243514.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

