

Product Datasheet

AP3B2 antibody (orb140169)

Description	Rabbit Polyclonal to AP3B2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	AP3B2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: AP3B2. Gene-IDs: 8120. Accessions: NP_004635.2, Q13367.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	119059
Uniprot ID	Q13367
NCBI	NP_001265441.1, NP_004635.2, NP_001265440.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

