

Product Datasheet

USP27X antibody (orb140173)

Description	Rabbit Polyclonal to USP27X.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	USP27X
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: USP27X. Gene-IDs: 389856. Accessions: A6NNY8, NP_001138545.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49630
Uniprot ID	A6NNY8
NCBI	NP_001138545.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

