

Product Datasheet

ABCF2 antibody (orb140175)

Description	Rabbit Polyclonal to ABCF2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ABCF2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ABCF2. Gene-IDs: 10061. Accessions: NP_009120.1, Q9UG63, NP_005683.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	71290
Uniprot ID	Q9UG63
NCBI	NP_005683.2, NP_009120.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

