

Product Datasheet

USP1 antibody (orb140176)

Description	Rabbit Polyclonal to USP1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	USP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: USP1. Gene-IDs: 7398. Accessions: NP_001017416.1, O94782, NP_001017415.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	88207
Uniprot ID	O94782
NCBI	NP_003359.3 , NP_001017415.1 , NP_001017416.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

