

Product Datasheet

UBA6 antibody (orb140178)

Description	Rabbit Polyclonal to UBA6.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	UBA6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: UBA6. Gene-IDs: 55236. Accessions: A0AVT1, NP_060697.4.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	117970
Uniprot ID	A0AVT1
NCBI	NP_060697.4
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

