

Product Datasheet

TBC1D8 antibody (orb140186)

Description	Rabbit Polyclonal to TBC1D8.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TBC1D8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: TBC1D8. Gene-IDs: 11138. Accessions: NP_001095896.1, O95759.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	130835
Uniprot ID	O95759
NCBI	NP_001095896.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

