

Product Datasheet

SEMA6D antibody (orb140187)

Description	Rabbit Polyclonal to SEMA6D.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SEMA6D
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: SEMA6D. Gene-IDs: 80031. Accessions: Q8NFY4, NP_065909.1, NP_079242.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	119872
Uniprot ID	Q8NFY4
NCBI	NP_001185928.1, NP_705870.1, NP_705869.1, NP_705871.1, NP_065909.1, NP_705872.1, NP_079242.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

