

Product Datasheet

KIF5B antibody (orb140190)

Description	Rabbit Polyclonal to KIF5B.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KIF5B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: KIF5B. Gene-IDs: 3799. Accessions: P33176, NP_004512.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	109685
Uniprot ID	P33176
NCBI	NP_004512.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

