

Product Datasheet

ChAT (Choline Acetyltransferase) Antibody (orb1401953)

Description
Our Anti-ChAT (Choline

from Biorbyt is goat polyclonal. It
detec...

Species/Host

Goat

Reactivity

Gallus, Guinea pig, Human, Mouse,
Primate, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

Recombinant human choline
acetyltransferase

Target

ChAT (Choline Acetyltransferase)

Preservatives

PBS + 5 mg/ml rAlbumin + 0.2% NaN₃

Form/Appearance

Affinity Purified

Storage

Recommended that the undiluted
antibody be aliquoted into smaller
working volumes (10-30 uL/vial
depending on usage) upon arrival and
stored long term at -20° C or -80° C,
while keeping a working aliquot stored
at 4° C for short term. Avoid
freeze/thaw cycles. Stable for at least
1 year.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

70

Uniprot ID

P28329

Dilution Range

WB: 1:100, IHC: 1:100-1:1000

Expiration Date

12 months from date of receipt.