

Product Datasheet

PSG2 antibody (orb140197)

Description	Rabbit Polyclonal to PSG2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PSG2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: PSG2. Gene-IDs: 5670. Accessions: NP_112536.2, P11465.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37216
Uniprot ID	P11465
NCBI	NP_112536.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

