

Product Datasheet

TSSC1 antibody (orb140199)

Description	Rabbit Polyclonal to TSSC1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	EIPR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: TSSC1. Gene-IDs: 7260. Accessions: Q53HC9, NP_003301.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43603
Uniprot ID	Q53HC9
NCBI	NP_003301.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

