

Product Datasheet

CSPG5 antibody (orb140200)

Description	Rabbit Polyclonal to CSPG5.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CSPG5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: CSPG5. Gene-IDs: 10675. Accessions: NP_006565.2, O95196.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60016
Uniprot ID	O95196
NCBI	NP_001193872.1 , NP_001193873.1 , NP_006565.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

