

Product Datasheet

TPRX1 antibody (orb140208)

Description	Rabbit Polyclonal to TPRX1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TPRX1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: TPRX1. Gene-IDs: 284355. Accessions: NP_940881.2, Q8N7U7.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	40603
Uniprot ID	Q8N7U7
NCBI	NP_940881.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

