

Product Datasheet

SELENOT Antibody (orb1402102)

Description

SELENOT Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC, WB

Applications
Immunogen

E.coli-derived human SELENOT recombinant protein
(Position: S20-S195).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of
500 µg/ml.

Storage

At -20°C for one year from date of receipt. After
reconstitution, at 4°C for one month. It can also be
aliquotted and stored frozen at -20°C for six months. Avoid
repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive
results. Other applications have not been tested. Optimal
dilutions should be determined by end users. Adding 0.2 ml
of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW




22 kDa

Uniprot ID
P62341
Dilution Range

Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat
Immunohistochemistry(Paraffin-embedded Section), 2-5
µg/ml, Human, Mouse, Rat Immunofluorescence, 5 µg/ml,
Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human Direct
ELISA, 0.1-0.5 µg/ml, Human

Expiration Date

12 months from date of receipt.


IHC analysis
of SELENOT
using anti-
SELEN...

IHC analysis
of SELENOT
using anti-
SELEN...

IHC analysis
of SELENOT
using anti-
SELEN...