

Product Datasheet

Lysozyme-like protein 6 antibody (orb140211)

Description	Rabbit polyclonal antibody to Lysozyme-like protein 6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LYZL6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: LYZL6. Gene-IDs: 57151. Accessions: O75951, NP_065159.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16956
Uniprot ID	O75951
NCBI	NP_065159.1 , NP_001186880.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

