

Product Datasheet

Selenophosphate synthetase 1/SEPHS1 Antibody (orb1402161)

Description

Selenophosphate synthetase 1/SEPHS1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence of human Selenophosphate synthetase 1/SEPHS1 (RERDKVMPLIIQGFKDAEEA).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

50 kDa

Uniprot ID



P49903

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat
Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Rat

Expiration Date

12 months from date of receipt.



IHC analysis of
Selenophosphate
Synthetase...



Western blot
analysis of
Selenophosphate...