

Product Datasheet

ATXN7 antibody (orb140218)

Description	Rabbit Polyclonal to ATXN7.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ATXN7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ATXN7. Gene-IDs: 6314. Accessions: NP_000324.1, O15265, NP_001170858.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	95451
Uniprot ID	O15265
NCBI	NP_000324.1, NP_001170858.1, NP_001121621.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

