

Product Datasheet

Seryl-tRNA synthetase/SERS/SARS1 Antibody (orb1402183)

Description

Seryl-tRNA synthetase/SERS/SARS1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, FC, ICC, IF, WB

Applications

Immunogen

E.coli-derived human Seryl-tRNA synthetase/SERS/SARS1 recombinant protein (Position: K71-K483).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

59 kDa

Uniprot ID



[P49591](#)

Dilution Range


Western blot, 0.1-0.25 µg/ml, Human, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 5 µg/ml,
Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human
Direct ELISA, 0.1-0.5 µg/ml, Human

Expiration Date


12 months from date of receipt.



IF analysis of
Seryl-TRNA
Synthetase/SER...



Western blot
analysis of Seryl-
TRNA Synt...



Flow Cytometry
analysis of Caco-
2 cells ...