






Product Datasheet

CD74 Antibody (orb1402197)

Description	CD74 Antibody	
Species/Host	Rabbit	 
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, FC, IF, IHC, WB	
Immunogen	E.coli-derived human CD74 recombinant protein (Position: Q73-M296).	IHC analysis of CD74 using anti-CD74 ant...
Form/Appearance	Lyophilized	 
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	IHC analysis of CD74 using anti-CD74 ant...
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	31-34 kDa	
Uniprot ID	P04233	
Dilution Range	Western blot, 0.25-0.5 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10 ⁶ cells, Human Direct ELISA, 0.1-0.5 µg/ml, Human	IF analysis of CD74 using anti-CD74 anti...
Expiration Date	12 months from date of receipt.	