


Product Datasheet

Human FGL1 ELISA Kit (orb1402205)

Description	Human FGL1 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L23-I312
Form/Appearance	Lyophilized powder
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q08830
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human FGL1
PicoKine
ELISA Kit
standard c...