

Product Datasheet

DEXI antibody (orb140230)

Description	Rabbit Polyclonal to DEXI.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DEXI
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: DEXI. Gene-IDs: 28955. Accessions: NP_054734.2, O95424.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	10429
Uniprot ID	O95424
NCBI	NP_054734.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

