

Product Datasheet

MRPL1 antibody (orb140236)

Description	Rabbit Polyclonal to MRPL1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MRPL1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: MRPL1. Gene-IDs: 65008. Accessions: Q9BYD6, NP_064621.3.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36909
Uniprot ID	Q9BYD6
NCBI	NP_064621.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

