

Product Datasheet

FKBP6 antibody (orb140241)

Description	Rabbit Polyclonal to FKBP6.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	FKBP6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: FKBP6. Gene-IDs: 8468. Accessions: O75344, NP_003593.3.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37214
Uniprot ID	O75344
NCBI	NP_003593.3, NP_001128683.1, NP_001268233.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

