

Product Datasheet

SRM antibody (orb140244)

Description	Rabbit Polyclonal to SRM.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SRM
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: SRM. Gene-IDs: 6723. Accessions: NP_003123.2, P19623.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33825
Uniprot ID	P19623
NCBI	NP_003123.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

