

Product Datasheet

MTIF2 antibody (orb140247)

Description	Rabbit Polyclonal to MTIF2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MTIF2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: MTIF2. Gene-IDs: 4528. Accessions: NP_001005369.1, P46199, NP_002444.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	81317
Uniprot ID	P46199
NCBI	NP_001005369.1, NP_002444.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

