

Product Datasheet

SRRD antibody (orb140252)

Description	Rabbit Polyclonal to SRRD.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SRRD
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: SRRD. Gene-IDs: 402055. Accessions: Q9UHH36, NP_001013716.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38573
Uniprot ID	Q9UHH36
NCBI	NP_001013716.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

