

Product Datasheet

IL1RAPL2 antibody (orb140255)

Description	Rabbit Polyclonal to IL1RAPL2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	IL1RAPL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: IL1RAPL2. Gene-IDs: 26280. Accessions: NP_059112.1, Q9NP60.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	78670
Uniprot ID	Q9NP60
NCBI	NP_059112.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

