

Product Datasheet

MAGEB16 antibody (orb140257)

Description	Rabbit Polyclonal to MAGEB16.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MAGEB16
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: MAGEB16. Gene-IDs: 139604. Accessions: A2A368, NP_001093391.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36178
Uniprot ID	A2A368
NCBI	NP_001093391.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

