

Product Datasheet

OR4K1 antibody (orb140258)

Description	Rabbit Polyclonal to OR4K1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	OR4K1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: OR4K1. Gene-IDs: 79544. Accessions: Q8NGD4, NP_001004063.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35201
Uniprot ID	Q8NGD4
NCBI	NP_001004063.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

