

Product Datasheet

RIMBP3 antibody (orb140260)

Description	Rabbit Polyclonal to RIMBP3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RIMBP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: RIMBP3. Gene-IDs: 85376. Accessions: Q9UFD9.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	180717
Uniprot ID	Q9UFD9
NCBI	NP_056487.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

