

Product Datasheet

KRT6A antibody (orb140265)

Description	Rabbit Polyclonal to KRT6A.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KRT6A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: KRT6A. Gene-IDs: 3853. Accessions: P02538, NP_005545.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60045
Uniprot ID	P02538
NCBI	NP_005545.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

