

Product Datasheet

PRAMEF14 antibody (orb140282)

Description	Rabbit Polyclonal to PRAMEF14.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PRAMEF14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: PRAMEF14. Gene-IDs: 729528. Accessions: Q5SWL7, NP_001093324.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55081
Uniprot ID	Q5SWL7
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

