

Product Datasheet

ACTL6A antibody (orb140292)

Description	Rabbit Polyclonal to ACTL6A.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ACTL6A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ACTL6A. Gene-IDs: 86. Accessions: NP_817126.1, O96019, NP_004292.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47461
Uniprot ID	O96019
NCBI	NP_829888.1, NP_817126.1, NP_004292.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

