

Product Datasheet

TNNT3 antibody (orb140309)

Description	Rabbit Polyclonal to TNNT3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TNNT3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: TNNT3. Gene-IDs: 7140. Accessions: NP_001036246.1, NP_001036245.1, P45378.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31825
Uniprot ID	P45378
NCBI	NP_001036246.1, NP_001036247.1, NP_006748.1, NP_001036245.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

