

Product Datasheet

RGS11 antibody (orb140312)

Description	Rabbit Polyclonal to RGS11.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RGS11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: RGS11. Gene-IDs: 8786. Accessions: NP_899180.1, O94810, NP_003825.1.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52946
Uniprot ID	O94810
NCBI	NP_899180.1, NP_003825.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

