

Product Datasheet

ZNF408 antibody (orb140328)

Description	Rabbit polyclonal antibody to ZNF408 also known as PRDM17 (PR domain zinc finger protein 17) has 720...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZNF408
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ZNF408. Gene-IDs: 79797. Accessions: NP_001171680.1, NP_079017.1, Q9H9D4.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	78439
Uniprot ID	Q9H9D4
NCBI	NP_001171680.1 , NP_079017.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

