

Product Datasheet

SNRNP70 antibody (orb140335)

Description	Rabbit Polyclonal to SNRNP70.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SNRNP70
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: SNRNP70. Gene-IDs: 6625. Accessions: NP_003080.2, P08621.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51557
Uniprot ID	P08621
NCBI	NP_003080.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

