

Product Datasheet

SAP30 antibody (orb140342)

Description	Rabbit Polyclonal to SAP30.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SAP30
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: SAP30. Gene-IDs: 8819. Accessions: NP_003855.1, O75446.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23306
Uniprot ID	O75446
NCBI	NP_003855.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

