

Product Datasheet

RRM2B antibody (orb140343)

Description	Rabbit Polyclonal to RRM2B.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RRM2B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: RRM2B. Gene-IDs: 50484. Accessions: NP_001165948.1, NP_001165949.1, Q7LG56.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	40737
Uniprot ID	Q7LG56
NCBI	NP_001165948.1, NP_001165949.1, NP_056528.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

