

Product Datasheet

CRLF1 antibody (orb140348)

Description	Rabbit Polyclonal to CRLF1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CRLF1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: CRLF1. Gene-IDs: 9244. Accessions: NP_004741.1, O75462.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46302
Uniprot ID	O75462
NCBI	NP_004741.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

