

Product Datasheet

MAP3K7IP3 antibody (orb140360)

Description	Rabbit Polyclonal to MAP3K7IP3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TAB3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: TAB3. Gene-IDs: 257397. Accessions: Q8N5C8, NP_690000.2.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	78653
Uniprot ID	Q8N5C8
NCBI	NP_690000.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

