

Product Datasheet

RAPGEF5 antibody (orb140362)

Description	Rabbit Polyclonal to RAPGEF5.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RAPGEF5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: RAPGEF5. Gene-IDs: 9771. Accessions: NP_036426.3, Q92565.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	67733
Uniprot ID	Q92565
NCBI	NP_036426.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

