

Product Datasheet

ZNF649 antibody (orb140365)

Description	Rabbit Polyclonal to ZNF649.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZNF649
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Application notes	Gene Symbol: ZNF649. Gene-IDs: 65251. Accessions: NP_075562.2, Q9BS31.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57683
Uniprot ID	Q9BS31
NCBI	NP_075562.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

